Mepilex® Border Flex Carré

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(2) Mölnlycke Health Care. Mepilex® Border Flex – External data – FEM simulations MxB Flex vs Allevyn Life. Report no. PD-529747. Data on file.

(3) Mölnlycke Health Care. FEM – Mepilex® Border Flex vs Aquacel foam. Report no. PD-534571. Data on file.

(4) Alten. Comparison of Mepilex® Border Flex dressing and Mepilex® Border dressing in wet condition. Report no. PD-530246. 18 January 2017. Data on file.

(5) ProDerm study report 16.0180-23. Assessment of Wearing Properties of Wound Dressings on the Knees. PD-535012.

(6) ProDerm study report 16.0180-23. Assessment of Wearing Properties of Wound Dressings on the Elbows. PD-535013.

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Mepilex® Border Flex

Self-adherent soft silicone foam dressing

Product description

Mepilex Border Flex is a self-adherent, absorbent dressing that maintains a moist wound environment. The waterproof outer layer protects the wound from dirt and bacteria. The dressing has a Safetac® wound contact layer that is a unique adhesive technology. It minimises pain to patients and trauma to wounds and the surrounding skin at dressing removal.

Mepilex Border Flex consists of:

- a wound contact layer consisting of soft silicone adhesive (Safetac) and a film carrier
- a flexible absorbent pad in three layers: a foam, a non-woven spreading layer and a layer with super absorbent fibres; the wound pad is partly perforated with Flex cut technology
- an outer film which is breathable but waterproof, providing a barrier to external contaminants

Dressing material content:

Silicone, polyurethane, polyacrylate, cotton, viscose, polyester and polyolefin.

Indications for use

Mepilex Border Flex is designed for a wide range of exuding wounds such as pressure ulcers, leg and foot ulcers, traumatic wounds (e.g. skin tears) and surgical wounds. Mepilex Border Flex can also be used on dry/necrotic wounds in combination with gels.

Mepilex Border Flex reduces postoperative blistering, and may also be used as part of a prophylactic therapy to help prevent skin damage, e.g. pressure ulcers.

Precautions

- Do not use on patients with known hypersensitivity to the ingoing materials/components of the product.
- Do not use together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- If you see signs of infection e.g. fever or the wound or surrounding skin becoming red, warm or swollen, consult a health care professional for appropriate treatment.
- The use of dressings as part of a prophylactic therapy does not preclude the need to continue to develop and follow a comprehensive pressure ulcer prevention protocol, i.e. support surfaces, positioning, nutrition, hydration, skin care and mobility.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if sterile barrier is damaged or opened prior to use. Do not re-sterilise.

Instructions for use

Mepilex Border Flex can be used by lay persons under supervision of health care professionals.

1. Select an appropriate dressing size/shape. When applicable, choose a dressing shape optimized for the specific anatomical site.

For use on wounds:

The wound pad should overlap the dry surrounding skin by at least 1-2cm in order to protect the surrounding skin from maceration and to affix the dressing securely.

If used for prevention:

Assure that the wound pad covers the area at risk of pressure injury/tissue damage.

2. For use on wounds:

Cleanse the wound in accordance to clinical practice. Dry the surrounding skin thoroughly.

3. Remove the first release film and apply the adherent side to the wound.

4. Remove the remaining release liners and smooth down the border on the skin.
Do not stretch.

Mepilex Border Flex may be left in place up to 7 days, depending on the condition of the wound and surrounding skin, or as indicated by clinical practice.

For use on wounds:

- Mepilex Border Flex can be used under compression bandaging.
- As Mepilex Border Flex maintains a moist wound environment, which supports debridement, there might be an initial increase in the wound size. This is normal and to be expected.
- A change in dressing regimen can result in an initial increased level of exudates, which temporarily may require an increased change frequency.

If used for prevention:

The area at risk of pressure ulcer/tissue damage should be inspected at regular intervals according to clinical practice.

Disposal should be handled according to local environmental procedures.

Other information

The polyurethane foam used in the product may change colour to more yellow when it is exposed to light, air and/or heat. The colour change has no influence on product properties when used before expiry date.

If any serious incident has occurred in relation to the use of Mepilex Border Flex it should be reported to Mölnlycke Health Care.

Mepilex® and Safetac® are registered trademarks of Mölnlycke Health Care AB.

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CANADA Mölnlycke Health Care Inc.,
2010 Winston Park Drive, Suite 100,
Oakville, Ontario, L6H 5R7
Tel: +1 905 829 1502

CHINA Mölnlycke Health Care
(Shanghai) Co., Ltd.
RM 629, No.8 Hua Jing Road,
Wai Gao Qiao FTZ Shanghai, China
Tel: +86 10 5128 8571

CZECH REPUBLIC Mölnlycke Health Care s.r.o.
Hajkova 2747/22, 130 00 Prague 3
Tel: +420 221 890 517 (**Reception**)
Tel: +420 221 890 511 (**Customer Service**)

DENMARK Mölnlycke Health Care ApS,
Gydevang 33, DK-3450 Allerød
Tel: +45 48 16 82 68
Customer Service Tel: +80 886 810

ESTONIA Mölnlycke Health Care OÜ,
Vabaõhumuuseumitee 4-67, EE-13522 Tallin
Tel: +372 671 1520

FINLAND Mölnlycke Health Care Oy,
Pitäjänmäentie 14, 3 krs, FI-00380 Helsinki
Tel: +358 201 622 300
Customer Service Tel: 0800 113 442

FRANCE Mölnlycke Health Care SAS,
13 Allée du Château Blanc
59290 WASQUEHAL
Tel: +33 320 122 555
Customer Service Tel: 0800 910 292

GERMANY Mölnlycke Health Care GmbH
Grafenberger Allee 297
DE-40237 Düsseldorf
Tel: +49 211 92 08 80
Customer Service Tel: 0800 186 21 80

HUNGARY Mölnlycke Health Care Kft
HU-1134 Budapest, Dévai u. 26-28. III. em.
Tel: +36 1 477 3080

INDIA Mölnlycke Health Care India Pvt Ltd
Universal Business Centre
Second Floor, Universal Trade Tower
Sector 49, Gurgaon Sohna Road
Gurgaon 122 018, Haryana
Tel: +91 124 4696 008

ITALY Mölnlycke Health Care s.r.l.,
Via Marsala 40/C, IT-21013 Gallarate (VA)
Tel: +39 0331 714 011

JAPAN Mölnlycke Health Care KK,
4F Conceria Nishi-Shinjuku Tower's West
6-20-7 Nishi-Shinjuku,
Shinjuku-Ku, Tokyo 160-0023
Tel: +81 3 6914 5004

KOREA Mölnlycke Health Care Korea Co., Ltd
2nd floor, Hansol Richville.
46 Wiryeseong-daero,
Songpa-gu, Seoul, 05627
South Korea
Tel: +82 2 3789 1402

LATVIA Mölnlycke Health Care AB,
Dzelzavas 120 M, Riga LV-1021
Tel: +372 6711 522

LITHUANIA Mölnlycke Health Care AB,
J. Galvydzio g. 3, LT 08236 Vilnius
Tel: +370 5 2274 58 37

NETHERLANDS Mölnlycke Health Care B.V.
Postbus 3196, 4800 DD Breda
Tel: +31 76 521 9663

NORWAY Mölnlycke Health Care AS
Postboks 6229 Etterstad, NO-0603 Oslo
Tel: +47 22 70 63 70
Customer Service Tel: 800 161 78

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USA 1-800-882-4582
Canada 1-800-494-5134

**EU Declaration of Conformity
according to MDR**

Document ID: PD-705293 **Rev:**1
Created by: Sofia Horkeby
State: Released
Approved by: Christina Lewing
Release date: 2021-11-01 09:15:51

Title: Mepilex Border Flex/Mepilex Border Comfort

Dates and times in Greenwich Mean Time, 24 hours format

We, Mölnlycke Health Care AB, Gamlestadsvägen 3C, Box 13080, SE-402 52 Göteborg, Sweden being exclusively responsible manufacturer for conformity of the device/s; declare that the devices listed in the attached schedule are in conformity with the provisions of the Medical Device Regulation 2017/745, concerning medical devices.

Other Union Legislation applicable: Not applicable

Trade name / Product name:	Mepilex Border Flex/Mepilex Border Comfort
----------------------------	--

Product Classification: MDR Class IIb

MDR Classification Rule: 4

Sterility: Sterile

Measuring Function: No

This declaration is supported by a conformity assessment procedure in accordance with Annex/es: IX

Common Specification: No CS is applicable

Certificate number: MDR 722028

Issued by: BSi

Id No: 2797

(Notified Body Name)

Note. For non-sterile, non-measuring class I products there is no Notified Body certificate

Mölnlycke Health Care issues this declaration in recognition of all applicable harmonised standards.

Signed for and on behalf of Mölnlycke Health Care

Place of Issue: Göteborg, Sweden

Date: 2021-11-01

Function: Global Regulatory Affairs Director

Name: Christina Lewing

Signature:

**EU Declaration of Conformity
according to MDR**

Document ID: PD-705293 **Rev:**1
Created by: Sofia Horkeby
State: Released
Approved by: Christina Lewing
Release date: 2021-11-01 09:15:51

Title: Mepilex Border Flex/Mepilex Border Comfort

Dates and times in Greenwich Mean Time, 24 hours format

Generic product information	
Intended Purpose:	Mepilex Border Flex is designed for a wide range of exuding wounds such as pressure ulcers, leg and foot ulcers, traumatic wounds (e.g. skin tears) and surgical wounds. Mepilex Border Flex can also be used on dry/necrotic wounds in combination with gels. Mepilex Border Flex reduces postoperative blistering, and may also be used as part of a prophylactic therapy to help prevent skin damage, e.g. pressure ulcers.
Basic UDI-DI:	733243000000000036JV
GMDN Code:	46854 Wound - nonadherent dressing, absorbent, sterile

Product References Covered by this Declaration:	Product Descriptor:
583300	Self-adherent soft silicone foam dressing
583320	Self-adherent soft silicone foam dressing
583400	Self-adherent soft silicone foam dressing
583420	Self-adherent soft silicone foam dressing
583500	Self-adherent soft silicone foam dressing
583520	Self-adherent soft silicone foam dressing
583620	Self-adherent soft silicone foam dressing
595000	Self-adherent soft silicone foam dressing
595011	Self-adherent soft silicone foam dressing
595021	Self-adherent soft silicone foam dressing
595200	Self-adherent soft silicone foam dressing
595211	Self-adherent soft silicone foam dressing
595221	Self-adherent soft silicone foam dressing
595277	Self-adherent soft silicone foam dressing
595300	Self-adherent soft silicone foam dressing
595311	Self-adherent soft silicone foam dressing
595321	Self-adherent soft silicone foam dressing
595377	Self-adherent soft silicone foam dressing

**EU Declaration of Conformity
according to MDR**

Document ID: PD-705293 **Rev:**1
Created by: Sofia Horkeby
State: Released
Approved by: Christina Lewing
Release date: 2021-11-01 09:15:51

Title: Mepilex Border Flex/Mepilex Border Comfort

Dates and times in Greenwich Mean Time, 24 hours format

595400	Self-adherent soft silicone foam dressing
595411	Self-adherent soft silicone foam dressing
595421	Self-adherent soft silicone foam dressing
595477	Self-adherent soft silicone foam dressing
595600	Self-adherent soft silicone foam dressing
595611	Self-adherent soft silicone foam dressing
595621	Self-adherent soft silicone foam dressing
595800	Self-adherent soft silicone foam dressing
595820	Self-adherent soft silicone foam dressing
595850	Self-adherent soft silicone foam dressing
595900	Self-adherent soft silicone foam dressing
595920	Self-adherent soft silicone foam dressing
596011	Self-adherent soft silicone foam dressing
596211	Self-adherent soft silicone foam dressing
596311	Self-adherent soft silicone foam dressing
596411	Self-adherent soft silicone foam dressing
596611	Self-adherent soft silicone foam dressing
596800	Self-adherent soft silicone foam dressing
596900	Self-adherent soft silicone foam dressing

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© Copyright (2020) Mölnlycke Health Care.
Mepilex® och Safetac® are registered trademarks of Mölnlycke Health Care AB.