

TECHNICAL DATA SHEET N° 13

074659P-P Product Description : *Medical Roll Microembossed glued*

Customer : **MAURY**

		Standard	Unit	Data	
1	Composition	Virgin pulp			
2	Grammage	UNI EN ISO 12625-6	g/m²	34	± 5%
3	Ply No.de plis	-	no.	2	
4	MD Tensile strength Rèst.a sec-Sens Marche	UNI EN ISO 12625-4	N/m	410	± 20%
5	CD Tensile strength Rèst.a sec-Sens Travers	UNI EN ISO 12625-4	N/m	240	± 20%
6	MD Wet tensile strength Rèst.mouillè-Sens Marche	UNI EN ISO 12625-5	N/m	82	± 20%
7	CD Wet tensile strength Rèst.mouillè-Sens Travers	UNI EN ISO 12625-5	N/m	48	± 20%
8	Brightness Gradation de blanc	UNI EN ISO 12625-7	%	85	± 3%
9	Roll Width Rouleau Hauteur	-	mm	590	± 5
10	Sheets Length Feuille Longueur	-	mm	340	± 5
11	Sheets count Formats	-	no.	135	± 3%
12	Roll length Longeur du rouleau	-	m	45,9	± 5%
13	Roll weight Poids	-	g	980	±5%
14	Inside core diameter Diamètre intérieur mandrin		mm	40	±1
15	Embossed Gaufrè	-	Microembossed		
16	Packing Emballage	pack	single PE bag / Sachet en PE		
		carton	9 rolls / rouleau		
		pallet	24 cartons / Carton		
Redatta da:		Verificata da:	Approvata da:	Revisione: R01	
Nome:C.LIBRIZZI		Nome:C.LIBRIZZI	Nome:I.RACCA		