

SPECIFICATIONS

ROMED STITCH CUTTER

REF: STCU-5000



Applicable harmonized standard: EN 27740:1992
Medical device conform MDD 93/42/EEC

Material: Carbon steel SK4 (Japanese standard, as Cr06)

Hardness: $\geq 750\text{HV}10$

Sharpness: $\leq 0.8\text{ N}$

Thickness: $0.40 \pm 0.02\text{ mm}$

Length: 6.5 cm

Sterilization: Gamma Irradiation

Expiry date: 5 years after sterilisation

Packing

Per piece in an aluminium foil, 100 pieces in a box, 50 boxes in an export carton